

# PERFADEX® Plus

The gold standard in lung preservation

PERFADEX® Plus is an extracellular, low potassium, dextran-based electrolyte preservation solution indicated for the flushing, cold static storage and transportation of isolated lungs after removal from the donor in preparation for eventual transplantation into a recipient.

- PERFADEX® Plus is **READY, EASY** and **SAFE** to use.
- PERFADEX® Plus **enables safe preservation of lungs up to 12 hours**, depending on the quality of the donor lung being retrieved.<sup>1</sup>
- **Gold Standard for lung preservation for >20 years**
- **Clinical evidence from more than 65.000 lung transplantations\***
- **Used by more than 90% of all lung transplantation clinics**

\*15,000 lung transplantations with PERFADEX® Plus and 50,000 with PERFADEX®

1. Yeung JC, et al. Outcomes after transplantation of lungs preserved for more than 12 h: a retrospective study. *Lancet Respir Med.* 2017 Feb;5(2):119-124.

## Order information:

### Reference no:

19811, 19850 (US)

19911, 19950 (US)

### Product:

10 x PERFADEX® Plus 1000 ml

10 x XVIVO Click Adapter

2 x PERFADEX® Plus 3000 ml

2 x XVIVO Click Adapter

XVIVO Perfusion AB, Mässans gata 10, SE-412 51 Göteborg, Sweden  
Tel +46 31 788 21 50 | [xvivogroup.com](http://xvivogroup.com)

PERFADEX® Plus is a trademark of XVIVO Perfusion AB.

© 2021 XVIVO Perfusion AB. All rights reserved.

Issued: 2021-10 Doc ID: MAR-16377-v.1.0



**XVIVO**